



The partner you trust

Bone Screws Catalogue

About Us

EgiFix Medical is a leading company in Egypt and the Middle East in the field of manufacturing orthopedic medical implants.

We are dedicated to providing surgeons with a wide variety of products with the aim of giving patients the ultimate care they deserve. EgiFix Medical offers innovative solutions in spine surgery, Arthroplasty and trauma fixation.

Cairo-based EgiFix was founded in 2003. With two manufacturing premises.

We have been writing a success story for 20 years - and counting.

INDEX

<u>Cortical Screws</u>	<u>4</u>
<u>Cancellous Screws</u>	<u>12</u>
<u>Cannulated Compression Screws</u>	<u>20</u>
<u>Cannulated Headless Compression Screws</u>	<u>28</u>
<u>Plate locking screws</u>	<u>32</u>

Cortical Screws

Design Features

- ISO Standard 5835
- Hexagonal head for strong and better alignment with screw driver
- Hemispherical countersink provides higher angulation through washer or plate hole maintaining concentric contact between the screw and the sides of the plate
- Fully threaded screw
- Thread shape according to ISO standard
- The screw is self tapping

Material

- SS 316LVM
- Tit Grade 5 Eli

Sizes

- 2.0 mm
 - Length: Starting from 6mm to 30 mm with 2mm increments
- 2.7 mm
 - Length Starting from 6mm to 40 mm with 2mm increments
- 3.5 mm
 - Length Starting from 10mm to 40 mm with 2mm increments
 - 40mm to 110 mm with 5mm increments
- 4.5 mm
 - Length Starting from 10mm to 80 mm with 2mm increments
 - 80mm to 120 mm with 5mm increments

Cortical Screw \varnothing 2.0 mm

Core	\varnothing 1.3 mm
Outer Diameter	\varnothing 2.0 mm
Head Diameter	\varnothing 4.0 mm
Head Length	1.8 mm
Drill Bit for thread Hole	\varnothing 1.5 mm
Tap	HA 2.0
Screw driver	Hexagonal
Hexagonal Size	1.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
5 mm	EGTS-01S-2.0-005	EGTS-01T-2.0-005
6 mm	EGTS-01S-2.0-006	EGTS-01T-2.0-006
8 mm	EGTS-01S-2.0-008	EGTS-01T-2.0-008
10 mm	EGTS-01S-2.0-010	EGTS-01T-2.0-010
12 mm	EGTS-01S-2.0-012	EGTS-01T-2.0-012
14 mm	EGTS-01S-2.0-014	EGTS-01T-2.0-014
16 mm	EGTS-01S-2.0-016	EGTS-01T-2.0-016
18 mm	EGTS-01S-2.0-018	EGTS-01T-2.0-018
20 mm	EGTS-01S-2.0-020	EGTS-01T-2.0-020
22 mm	EGTS-01S-2.0-022	EGTS-01T-2.0-022
24 mm	EGTS-01S-2.0-024	EGTS-01T-2.0-024
26 mm	EGTS-01S-2.0-026	EGTS-01T-2.0-026
28 mm	EGTS-01S-2.0-028	EGTS-01T-2.0-028
30 mm	EGTS-01S-2.0-030	EGTS-01T-2.0-030

Cortical Screw \varnothing 2.7 mm

Core	\varnothing 1.9 mm
Outer Diameter	\varnothing 2.7 mm
Head Diameter	\varnothing 5.0 mm
Head Length	2.1 mm
Drill Bit for thread Hole	\varnothing 2.0 mm
Tap	HA 2.7
Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
6 mm	EGTS-01S-2.7-006	EGTS-01T-2.7-006
8 mm	EGTS-01S-2.7-008	EGTS-01T-2.7-008
10 mm	EGTS-01S-2.7-010	EGTS-01T-2.7-010
12 mm	EGTS-01S-2.7-012	EGTS-01T-2.7-012
14 mm	EGTS-01S-2.7-014	EGTS-01T-2.7-014
16 mm	EGTS-01S-2.7-016	EGTS-01T-2.7-016
18 mm	EGTS-01S-2.7-018	EGTS-01T-2.7-018
20 mm	EGTS-01S-2.7-020	EGTS-01T-2.7-020
22 mm	EGTS-01S-2.7-022	EGTS-01T-2.7-022
24 mm	EGTS-01S-2.7-024	EGTS-01T-2.7-024
26 mm	EGTS-01S-2.7-026	EGTS-01T-2.7-026
28 mm	EGTS-01S-2.7-028	EGTS-01T-2.7-028
30 mm	EGTS-01S-2.7-030	EGTS-01T-2.7-030
32 mm	EGTS-01S-2.7-032	EGTS-01T-2.7-032
34 mm	EGTS-01S-2.7-034	EGTS-01T-2.7-034
36 mm	EGTS-01S-2.7-036	EGTS-01T-2.7-036
38 mm	EGTS-01S-2.7-038	EGTS-01T-2.7-038
40 mm	EGTS-01S-2.7-040	EGTS-01T-2.7-040

Cortical Screw \varnothing 3.5 mm



Core	\varnothing 2.4 mm
Outer Diameter	\varnothing 3.5 mm
Head Diameter	\varnothing 6.0 mm
Head Length	3.1 mm
Drill Bit for thread Hole	\varnothing 2.5 mm
Tap	HA 3.5
Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
10 mm	EGTS-01S-3.5-010	EGTS-01T-3.5-010
12 mm	EGTS-01S-3.5-012	EGTS-01T-3.5-012
14 mm	EGTS-01S-3.5-014	EGTS-01T-3.5-014
16 mm	EGTS-01S-3.5-016	EGTS-01T-3.5-016
18 mm	EGTS-01S-3.5-018	EGTS-01T-3.5-018
20 mm	EGTS-01S-3.5-020	EGTS-01T-3.5-020
22 mm	EGTS-01S-3.5-022	EGTS-01T-3.5-022
24 mm	EGTS-01S-3.5-024	EGTS-01T-3.5-024
26 mm	EGTS-01S-3.5-026	EGTS-01T-3.5-026
28 mm	EGTS-01S-3.5-028	EGTS-01T-3.5-028
30 mm	EGTS-01S-3.5-030	EGTS-01T-3.5-030
32 mm	EGTS-01S-3.5-032	EGTS-01T-3.5-032
34 mm	EGTS-01S-3.5-034	EGTS-01T-3.5-034
36 mm	EGTS-01S-3.5-036	EGTS-01T-3.5-036
38 mm	EGTS-01S-3.5-038	EGTS-01T-3.5-038
40 mm	EGTS-01S-3.5-040	EGTS-01T-3.5-040
45 mm	EGTS-01S-3.5-045	EGTS-01T-3.5-045
50 mm	EGTS-01S-3.5-050	EGTS-01T-3.5-050
55 mm	EGTS-01S-3.5-055	EGTS-01T-3.5-055
60 mm	EGTS-01S-3.5-060	EGTS-01T-3.5-060
65 mm	EGTS-01S-3.5-065	EGTS-01T-3.5-065

Screw Length	S.S. Ref.	Tit. Ref
70 mm	EGTS-01S-3.5-070	EGTS-01T-3.5-070
75 mm	EGTS-01S-3.5-075	EGTS-01T-3.5-075
80 mm	EGTS-01S-3.5-080	EGTS-01T-3.5-080
85 mm	EGTS-01S-3.5-085	EGTS-01T-3.5-085
90 mm	EGTS-01S-3.5-090	EGTS-01T-3.5-090
95 mm	EGTS-01S-3.5-095	EGTS-01T-3.5-095
100 mm	EGTS-01S-3.5-100	EGTS-01T-3.5-100
105 mm	EGTS-01S-3.5-105	EGTS-01T-3.5-105
110 mm	EGTS-01S-3.5-110	EGTS-01T-3.5-110

Cortical Screw \varnothing 4.5 mm



Core	\varnothing 3.0 mm
Outer Diameter	\varnothing 4.5 mm
Head Diameter	\varnothing 8.0 mm
Head Length	4.6 mm
Drill Bit for thread Hole	\varnothing 3.2 mm
Tap	HA 4.5
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
10 mm	EGTS-01S-4.5-010	EGTS-01T-4.5-010
12 mm	EGTS-01S-4.5-012	EGTS-01T-4.5-012
14 mm	EGTS-01S-4.5-014	EGTS-01T-4.5-014
16 mm	EGTS-01S-4.5-016	EGTS-01T-4.5-016
18 mm	EGTS-01S-4.5-018	EGTS-01T-4.5-018
20 mm	EGTS-01S-4.5-020	EGTS-01T-4.5-020
22 mm	EGTS-01S-4.5-022	EGTS-01T-4.5-022
24 mm	EGTS-01S-4.5-024	EGTS-01T-4.5-024
26 mm	EGTS-01S-4.5-026	EGTS-01T-4.5-026
28 mm	EGTS-01S-4.5-028	EGTS-01T-4.5-028
30 mm	EGTS-01S-4.5-030	EGTS-01T-4.5-030
32 mm	EGTS-01S-4.5-032	EGTS-01T-4.5-032
34 mm	EGTS-01S-4.5-034	EGTS-01T-4.5-034
36 mm	EGTS-01S-4.5-036	EGTS-01T-4.5-036
38 mm	EGTS-01S-4.5-038	EGTS-01T-4.5-038
40 mm	EGTS-01S-4.5-040	EGTS-01T-4.5-040
42 mm	EGTS-01S-4.5-042	EGTS-01T-4.5-042
44 mm	EGTS-01S-4.5-044	EGTS-01T-4.5-044
46 mm	EGTS-01S-4.5-046	EGTS-01T-4.5-046
48 mm	EGTS-01S-4.5-048	EGTS-01T-4.5-048
50 mm	EGTS-01S-4.5-050	EGTS-01T-4.5-050
52 mm	EGTS-01S-4.5-052	EGTS-01T-4.5-052
54 mm	EGTS-01S-4.5-054	EGTS-01T-4.5-054
56 mm	EGTS-01S-4.5-056	EGTS-01T-4.5-056

Screw Length	S.S. Ref.	Tit. Ref
58 mm	EGTS-01S-4.5-058	EGTS-01T-4.5-058
60 mm	EGTS-01S-4.5-060	EGTS-01T-4.5-060
62 mm	EGTS-01S-4.5-062	EGTS-01T-4.5-062
64 mm	EGTS-01S-4.5-064	EGTS-01T-4.5-064
66 mm	EGTS-01S-4.5-066	EGTS-01T-4.5-066
68 mm	EGTS-01S-4.5-068	EGTS-01T-4.5-068
70 mm	EGTS-01S-4.5-070	EGTS-01T-4.5-070
72 mm	EGTS-01S-4.5-072	EGTS-01T-4.5-072
74 mm	EGTS-01S-4.5-074	EGTS-01T-4.5-074
76 mm	EGTS-01S-4.5-076	EGTS-01T-4.5-076
78 mm	EGTS-01S-4.5-078	EGTS-01T-4.5-078
80 mm	EGTS-01S-4.5-080	EGTS-01T-4.5-080
85 mm	EGTS-01S-4.5-085	EGTS-01T-4.5-085
90 mm	EGTS-01S-4.5-090	EGTS-01T-4.5-090
95 mm	EGTS-01S-4.5-095	EGTS-01T-4.5-095
100 mm	EGTS-01S-4.5-100	EGTS-01T-4.5-100
105 mm	EGTS-01S-4.5-105	EGTS-01T-4.5-105
110 mm	EGTS-01S-4.5-110	EGTS-01T-4.5-110
115 mm	EGTS-01S-4.5-115	EGTS-01T-4.5-115
120 mm	EGTS-01S-4.5-120	EGTS-01T-4.5-120



Cancellous Screws

Design Features

- ISO Standard 5835
- Hexagonal head for strong and better alignment with screw driver
- Hemispherical countersink provides higher angulation through washer or plate hole maintaining concentric contact between the screw and the sides of the plate
- Fully threaded screw
- Thread shape according to ISO standard
- The screw is self tapping

Material

- SS 316LVM
- Tit Grade 5 Eli

Sizes

- Diameter 4.0mm full thread
 - Length Starting from 12mm to 30 mm with 2mm increments
30mm to 70 mm with 5mm increments
- Diameter 4.0mm partial thread
 - Length Starting from 18mm to 40 mm with 2mm increments
40mm to 70 mm with 5mm increments
- Diameter 6.5mm full thread
 - Length Starting from 25mm to 120 mm with 5mm increments
- Diameter 6.5mm 16mm and 32 mm thread
 - Length Starting from 25mm to 120 mm with 5mm increments

Cancellous Screw \varnothing 4.0 Fully Threaded



Core	\varnothing 1.8 mm
Outer Diameter	\varnothing 4.0 mm
Head Diameter	\varnothing 6.0 mm
Head Length	2.7 mm
Drill Bit for thread Hole	\varnothing 2.0 mm
Tap	HB 4.0
Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
12 mm	EGTS-02S-4.0-012	EGTS-02T-4.0-012
14 mm	EGTS-02S-4.0-014	EGTS-02T-4.0-014
16 mm	EGTS-02S-4.0-016	EGTS-02T-4.0-016
18 mm	EGTS-02S-4.0-018	EGTS-02T-4.0-018
20 mm	EGTS-02S-4.0-020	EGTS-02T-4.0-020
22 mm	EGTS-02S-4.0-022	EGTS-02T-4.0-022
24 mm	EGTS-02S-4.0-024	EGTS-02T-4.0-024
26 mm	EGTS-02S-4.0-026	EGTS-02T-4.0-026
28 mm	EGTS-02S-4.0-028	EGTS-02T-4.0-028
30 mm	EGTS-02S-4.0-030	EGTS-02T-4.0-030
35 mm	EGTS-02S-4.0-035	EGTS-02T-4.0-035
40 mm	EGTS-02S-4.0-040	EGTS-02T-4.0-040
45 mm	EGTS-02S-4.0-045	EGTS-02T-4.0-045
50 mm	EGTS-02S-4.0-050	EGTS-02T-4.0-050
55 mm	EGTS-02S-4.0-055	EGTS-02T-4.0-055
60 mm	EGTS-02S-4.0-060	EGTS-02T-4.0-060

Cancellous Screw \varnothing 4.0 Partially Threaded



Core	\varnothing 1.8 mm
Outer Diameter	\varnothing 4.0 mm
Head Diameter	\varnothing 6.0 mm
Head Length	3.1 mm
Drill Bit for thread Hole	\varnothing 2.0 mm
Tap	HB 4.0
Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
18 mm	7 mm	EGTS-03S-4.0-018	EGTS-03T-4.0-018
20 mm	8 mm	EGTS-03S-4.0-020	EGTS-03T-4.0-020
22 mm	9 mm	EGTS-03S-4.0-022	EGTS-03T-4.0-022
24 mm	10 mm	EGTS-03S-4.0-024	EGTS-03T-4.0-024
26 mm	12 mm	EGTS-03S-4.0-026	EGTS-03T-4.0-026
28 mm	14 mm	EGTS-03S-4.0-028	EGTS-03T-4.0-028
30 mm	14 mm	EGTS-03S-4.0-030	EGTS-03T-4.0-030
35 mm	14 mm	EGTS-03S-4.0-035	EGTS-03T-4.0-035
40 mm	14 mm	EGTS-03S-4.0-040	EGTS-03T-4.0-040
45 mm	15 mm	EGTS-03S-4.0-045	EGTS-03T-4.0-045
50 mm	16 mm	EGTS-03S-4.0-050	EGTS-03T-4.0-050
55 mm	17 mm	EGTS-03S-4.0-055	EGTS-03T-4.0-055
60 mm	19 mm	EGTS-03S-4.0-060	EGTS-03T-4.0-060
65 mm	20 mm	EGTS-03S-4.0-065	EGTS-03T-4.0-065
70 mm	22 mm	EGTS-03S-4.0-070	EGTS-03T-4.0-070

Cancellous Screw \varnothing 6.5 Fully Threaded



Core	\varnothing 3.0 m
Outer Diameter	\varnothing 6.5 mm
Head Diameter	\varnothing 8.0 mm
Head Length	\varnothing 4.6 mm
Drill Bit for thread Hole	\varnothing 3.2 mm
Tap	HB 6.5
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Screw Length	S.S. Ref.	Tit. Ref
25 mm	EGTS-02S-6.5-025	EGTS-02T-6.5-025
30 mm	EGTS-02S-6.5-030	EGTS-02T-6.5-030
35 mm	EGTS-02S-6.5-035	EGTS-02T-6.5-035
40 mm	EGTS-02S-6.5-040	EGTS-02T-6.5-040
45 mm	EGTS-02S-6.5-045	EGTS-02T-6.5-045
50 mm	EGTS-02S-6.5-050	EGTS-02T-6.5-050
55 mm	EGTS-02S-6.5-055	EGTS-02T-6.5-055
60 mm	EGTS-02S-6.5-060	EGTS-02T-6.5-060
65 mm	EGTS-02S-6.5-065	EGTS-02T-6.5-065
70 mm	EGTS-02S-6.5-070	EGTS-02T-6.5-070
75 mm	EGTS-02S-6.5-075	EGTS-02T-6.5-075
80 mm	EGTS-02S-6.5-080	EGTS-02T-6.5-080
85 mm	EGTS-02S-6.5-085	EGTS-02T-6.5-085
90 mm	EGTS-02S-6.5-090	EGTS-02T-6.5-090
95 mm	EGTS-02S-6.5-095	EGTS-02T-6.5-095
100 mm	EGTS-02S-6.5-100	EGTS-02T-6.5-100
105 mm	EGTS-02S-6.5-105	EGTS-02T-6.5-105
110 mm	EGTS-02S-6.5-110	EGTS-02T-6.5-110
115 mm	EGTS-02S-6.5-115	EGTS-02T-6.5-115
120 mm	EGTS-02S-6.5-120	EGTS-02T-6.5-120

Cancellous Screw \varnothing 6.5 mm - 16 mm Threaded



Core	\varnothing 3.0 mm
Outer Diameter	\varnothing 6.5 mm
Head Diameter	\varnothing 8.0 mm
Head Length	\varnothing 4.6 mm
Shaft	\varnothing 4.5 mm
Drill Bit for thread Hole	\varnothing 3.2 mm
Tap	HB 6.5
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
25 mm	16 mm	EGTS-05S-6.5-025	EGTS-05T-6.5-025
30 mm	16 mm	EGTS-05S-6.5-030	EGTS-05T-6.5-030
35 mm	16 mm	EGTS-05S-6.5-035	EGTS-05T-6.5-035
40 mm	16 mm	EGTS-05S-6.5-040	EGTS-05T-6.5-040
45 mm	16 mm	EGTS-05S-6.5-045	EGTS-05T-6.5-045
50 mm	16 mm	EGTS-05S-6.5-050	EGTS-05T-6.5-050
55 mm	16 mm	EGTS-05S-6.5-055	EGTS-05T-6.5-055
60 mm	16 mm	EGTS-05S-6.5-060	EGTS-05T-6.5-060
65 mm	16 mm	EGTS-05S-6.5-065	EGTS-05T-6.5-065
70 mm	16 mm	EGTS-05S-6.5-070	EGTS-05T-6.5-070
75 mm	16 mm	EGTS-05S-6.5-075	EGTS-05T-6.5-075
80 mm	16 mm	EGTS-05S-6.5-080	EGTS-05T-6.5-080
85 mm	16 mm	EGTS-05S-6.5-085	EGTS-05T-6.5-085
90 mm	16 mm	EGTS-05S-6.5-090	EGTS-05T-6.5-090
95 mm	16 mm	EGTS-05S-6.5-095	EGTS-05T-6.5-095
100 mm	16 mm	EGTS-05S-6.5-100	EGTS-05T-6.5-100
105 mm	16 mm	EGTS-05S-6.5-105	EGTS-05T-6.5-105
110 mm	16 mm	EGTS-05S-6.5-110	EGTS-05T-6.5-110
115 mm	16 mm	EGTS-05S-6.5-115	EGTS-05T-6.5-115
120 mm	16 mm	EGTS-05S-6.5-120	EGTS-05T-6.5-120

Cancellous Screw \varnothing 6.5 mm - 32 mm Thread



Core	\varnothing 3.0 mm
Outer Diameter	\varnothing 6.5 mm
Head Diameter	\varnothing 8.0 mm
Head Length	\varnothing 4.6 mm
Shaft	\varnothing 4.5 mm
Drill Bit for thread Hole	\varnothing 3.2 mm
Tap	HB 6.5
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
35 mm	32 mm	EGTS-06S-6.5-035	EGTS-06T-6.5-035
40 mm	32 mm	EGTS-06S-6.5-040	EGTS-06T-6.5-040
45 mm	32 mm	EGTS-06S-6.5-045	EGTS-06T-6.5-045
50 mm	32 mm	EGTS-06S-6.5-050	EGTS-06T-6.5-050
55 mm	32 mm	EGTS-06S-6.5-055	EGTS-06T-6.5-055
60 mm	32 mm	EGTS-06S-6.5-060	EGTS-06T-6.5-060
65 mm	32 mm	EGTS-06S-6.5-065	EGTS-06T-6.5-065
70 mm	32 mm	EGTS-06S-6.5-070	EGTS-06T-6.5-070
75 mm	32 mm	EGTS-06S-6.5-075	EGTS-06T-6.5-075
80 mm	32 mm	EGTS-06S-6.5-080	EGTS-06T-6.5-080
85 mm	32 mm	EGTS-06S-6.5-085	EGTS-06T-6.5-085
90 mm	32 mm	EGTS-06S-6.5-090	EGTS-06T-6.5-090
95 mm	32 mm	EGTS-06S-6.5-095	EGTS-06T-6.5-095
100 mm	32 mm	EGTS-06S-6.5-100	EGTS-06T-6.5-100
105 mm	32 mm	EGTS-06S-6.5-105	EGTS-06T-6.5-105
110 mm	32 mm	EGTS-06S-6.5-110	EGTS-06T-6.5-110
115 mm	32 mm	EGTS-06S-6.5-115	EGTS-06T-6.5-115
120 mm	32 mm	EGTS-06S-6.5-120	EGTS-06T-6.5-120

Malleolar Cancellous ø 4.5 - Half Thread



Core	ø 3.0 mm
Outer Diameter	ø 4.5 mm
Head Diameter	ø 8.0 mm
Head Length	5.0 mm
Drill Bit for thread Hole	ø 3.2 mm
Tap	HA 4.5
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping - Trocar

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
20 mm	10 mm	EGTS-07S-4.5-020	EGTS-07T-4.5-020
25 mm	12 mm	EGTS-07S-4.5-025	EGTS-07T-4.5-025
30 mm	15 mm	EGTS-07S-4.5-030	EGTS-07T-4.5-030
35 mm	17 mm	EGTS-07S-4.5-035	EGTS-07T-4.5-035
40 mm	20 mm	EGTS-07S-4.5-040	EGTS-07T-4.5-040
45 mm	22 mm	EGTS-07S-4.5-045	EGTS-07T-4.5-045
50 mm	25 mm	EGTS-07S-4.5-050	EGTS-07T-4.5-050
55 mm	27 mm	EGTS-07S-4.5-055	EGTS-07T-4.5-055
60 mm	30 mm	EGTS-07S-4.5-060	EGTS-07T-4.5-060
65 mm	32 mm	EGTS-07S-4.5-065	EGTS-07T-4.5-065
70 mm	35 mm	EGTS-07S-4.5-070	EGTS-07T-4.5-070



Cannulated Screws

Design Features

- ISO Standard 5835
- Cannulated screw for insertion on guidewire
- Hexagonal head for strong and better alignment with screw driver
- Hemispherical countersink provides higher angulation through washer or plate hole maintaining concentric contact between the screw and the sides of the plate
- Partial threaded screw for better compression
- Thread shape according to ISO standard
- The screw is self drilling tip
- Cannulated instrumentation for easy use
- Reverse cutting distal thread

Material

- SS 316LVM
- Tit Grade 5 Eli

Sizes

- Diameter 4.0mm Short Thread
 - Length Starting from 25mm to 75 mm with 5mm increments
- Diameter 4.5mm Short thread
 - Length Starting from 10mm to 30 mm with 2mm increments
 - 30mm to 80 mm with 5mm increments
- Diameter 6.5mm (16/32 mm thread)
 - Length Starting from 40mm to 130 mm with 5mm increments
- Diameter 7.3mm (16/32 mm thread)
 - Length Starting from 40mm to 130 mm with 5mm increments

Cannulated Cancellous \varnothing 4.0 Partially Threaded



Core	\varnothing 2.4 mm
Outer Diameter	\varnothing 4.0 mm
Head Diameter	\varnothing 6.0 mm
Head Length	3.1 mm
Drill Bit for thread Hole	2.5 / 1.2
Guide Wire	\varnothing 1.0 mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
10 mm	4 mm	EGTS-09S-4.0-010	EGTS-09T-4.0-010
12 mm	4 mm	EGTS-09S-4.0-012	EGTS-09T-4.0-012
14 mm	5 mm	EGTS-09S-4.0-014	EGTS-09T-4.0-014
16 mm	5 mm	EGTS-09S-4.0-016	EGTS-09T-4.0-016
18 mm	6 mm	EGTS-09S-4.0-018	EGTS-09T-4.0-018
20 mm	7 mm	EGTS-09S-4.0-020	EGTS-09T-4.0-020
22 mm	7 mm	EGTS-09S-4.0-022	EGTS-09T-4.0-022
24 mm	8 mm	EGTS-09S-4.0-024	EGTS-09T-4.0-024
26 mm	8 mm	EGTS-09S-4.0-026	EGTS-09T-4.0-026
28 mm	9 mm	EGTS-09S-4.0-028	EGTS-09T-4.0-028
30 mm	10 mm	EGTS-09S-4.0-030	EGTS-09T-4.0-030
32 mm	11 mm	EGTS-09S-4.0-032	EGTS-09T-4.0-032
34 mm	11 mm	EGTS-09S-4.0-034	EGTS-09T-4.0-034
36 mm	12 mm	EGTS-09S-4.0-036	EGTS-09T-4.0-036
38 mm	12 mm	EGTS-09S-4.0-038	EGTS-09T-4.0-038
40 mm	13 mm	EGTS-09S-4.0-040	EGTS-09T-4.0-040
42 mm	14 mm	EGTS-09S-4.0-042	EGTS-09T-4.0-042
44 mm	14 mm	EGTS-09S-4.0-044	EGTS-09T-4.0-044
46 mm	15 mm	EGTS-09S-4.0-046	EGTS-09T-4.0-046
48 mm	15 mm	EGTS-09S-4.0-048	EGTS-09T-4.0-048
50 mm	16 mm	EGTS-09S-4.0-050	EGTS-09T-4.0-050

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
52 mm	16 mm	EGTS-09S-4.0-052	EGTS-09T-4.0-052
54 mm	17 mm	EGTS-09S-4.0-054	EGTS-09T-4.0-054
56 mm	18 mm	EGTS-09S-4.0-056	EGTS-09T-4.0-056
58 mm	19 mm	EGTS-09S-4.0-058	EGTS-09T-4.0-058
60 mm	20 mm	EGTS-09S-4.0-060	EGTS-09T-4.0-060
65 mm	22 mm	EGTS-09S-4.0-065	EGTS-09T-4.0-065
70 mm	24 mm	EGTS-09S-4.0-070	EGTS-09T-4.0-070
75 mm	26mm	EGTS-09S-4.0-075	EGTS-09T-4.0-075
80 mm	28 mm	EGTS-09S-4.0-080	EGTS-09T-4.0-080
85 mm	30 mm	EGTS-09S-4.0-085	EGTS-09T-4.0-085
90 mm	32 mm	EGTS-09S-4.0-090	EGTS-09T-4.0-090

Cannulated Cancellous \varnothing 4.5 Partially Threaded



Core	\varnothing 3.0 mm
Outer Diameter	\varnothing 4.5 mm
Head Diameter	\varnothing 6.0 mm
Head Length	4.6 mm
Drill Bit for thread Hole	3.2/1.2
Guide wire	\varnothing 1.6 mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping / Self Drilling

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
10 mm	5 mm	EGTS-09S-4.5-010	EGTS-09T-4.5-010
12 mm	5 mm	EGTS-09S-4.5-012	EGTS-09T-4.5-012
14 mm	6 mm	EGTS-09S-4.5-014	EGTS-09T-4.5-014
16 mm	12 mm	EGTS-09S-4.5-016	EGTS-09T-4.5-016
18 mm	12 mm	EGTS-09S-4.5-018	EGTS-09T-4.5-018
20 mm	12 mm	EGTS-09S-4.5-020	EGTS-09T-4.5-020
22 mm	12 mm	EGTS-09S-4.5-022	EGTS-09T-4.5-022
24 mm	12 mm	EGTS-09S-4.5-024	EGTS-09T-4.5-024
26 mm	12 mm	EGTS-09S-4.5-026	EGTS-09T-4.5-026
28 mm	12 mm	EGTS-09S-4.5-028	EGTS-09T-4.5-028
30 mm	12 mm	EGTS-09S-4.5-030	EGTS-09T-4.5-030
35 mm	12 mm	EGTS-09S-4.5-035	EGTS-09T-4.5-035
40 mm	12 mm	EGTS-09S-4.5-040	EGTS-09T-4.5-040
45 mm	12 mm	EGTS-09S-4.5-045	EGTS-09T-4.5-045
50 mm	16 mm	EGTS-09S-4.5-050	EGTS-09T-4.5-050
55 mm	16 mm	EGTS-09S-4.5-055	EGTS-09T-4.5-055
60 mm	16 mm	EGTS-09S-4.5-060	EGTS-09T-4.5-060
65 mm	16 mm	EGTS-09S-4.5-065	EGTS-09T-4.5-065
70 mm	16 mm	EGTS-09S-4.5-070	EGTS-09T-4.5-070
75 mm	16 mm	EGTS-09S-4.5-075	EGTS-09T-4.5-075
80 mm	16 mm	EGTS-09S-4.5-080	EGTS-09T-4.5-080

Canulated Cancellous $\phi 6.5$ -16mm thread



Core	$\phi 3.0$ mm
Outer Diameter	$\phi 6.5$ mm
Head Diameter	$\phi 8.0$ mm
Head Length	$\phi 4.6$ mm
Drill Bit for thread Hole	4.5/2.1
Guide wire	$\phi 2.0$ mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping / Self Drilling

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
40 mm	16 mm	EGTS-10S-6.5-040	EGTS-10T-6.5-040
45 mm	16 mm	EGTS-10S-6.5-045	EGTS-10T-6.5-045
50 mm	16 mm	EGTS-10S-6.5-050	EGTS-10T-6.5-050
55 mm	16 mm	EGTS-10S-6.5-055	EGTS-10T-6.5-055
60 mm	16 mm	EGTS-10S-6.5-060	EGTS-10T-6.5-060
65 mm	16 mm	EGTS-10S-6.5-065	EGTS-10T-6.5-065
70 mm	16 mm	EGTS-10S-6.5-070	EGTS-10T-6.5-070
75 mm	16 mm	EGTS-10S-6.5-075	EGTS-10T-6.5-075
80 mm	16 mm	EGTS-10S-6.5-080	EGTS-10T-6.5-080
85 mm	16 mm	EGTS-10S-6.5-085	EGTS-10T-6.5-085
90 mm	16 mm	EGTS-10S-6.5-090	EGTS-10T-6.5-090
95 mm	16 mm	EGTS-10S-6.5-095	EGTS-10T-6.5-095
100 mm	16 mm	EGTS-10S-6.5-100	EGTS-10T-6.5-100
105 mm	16 mm	EGTS-10S-6.5-105	EGTS-10T-6.5-105
110 mm	16 mm	EGTS-10S-6.5-110	EGTS-10T-6.5-110
115 mm	16 mm	EGTS-10S-6.5-115	EGTS-10T-6.5-115
120 mm	16 mm	EGTS-10S-6.5-120	EGTS-10T-6.5-120
125 mm	16 mm	EGTS-10S-6.5-125	EGTS-10T-6.5-125
130 mm	16 mm	EGTS-10S-6.5-130	EGTS-10T-6.5-130

Cannulated Cancellous $\phi 6.5$ -32mm thread



Core	$\phi 3.0$ mm
Outer Diameter	$\phi 6.5$ mm
Head Diameter	$\phi 8.0$ mm
Head Length	$\phi 4.6$ mm
Drill Bit for thread Hole	4.5/2.1
Guide wire	$\phi 2.0$ mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping / Self Drilling

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
40 mm	32 mm	EGTS-11S-6.5-040	EGTS-11T-6.5-040
45 mm	32 mm	EGTS-11S-6.5-045	EGTS-11T-6.5-045
50 mm	32 mm	EGTS-11S-6.5-050	EGTS-11T-6.5-050
55 mm	32 mm	EGTS-11S-6.5-055	EGTS-11T-6.5-055
60 mm	32 mm	EGTS-11S-6.5-060	EGTS-11T-6.5-060
65 mm	32 mm	EGTS-11S-6.5-065	EGTS-11T-6.5-065
70 mm	32 mm	EGTS-11S-6.5-070	EGTS-11T-6.5-070
75 mm	32 mm	EGTS-11S-6.5-075	EGTS-11T-6.5-075
80 mm	32 mm	EGTS-11S-6.5-080	EGTS-11T-6.5-080
85 mm	32 mm	EGTS-11S-6.5-085	EGTS-11T-6.5-085
90 mm	32 mm	EGTS-11S-6.5-090	EGTS-11T-6.5-090
95 mm	32 mm	EGTS-11S-6.5-095	EGTS-11T-6.5-095
100 mm	32 mm	EGTS-11S-6.5-100	EGTS-11T-6.5-100
105 mm	32 mm	EGTS-11S-6.5-105	EGTS-11T-6.5-105
110 mm	32 mm	EGTS-11S-6.5-110	EGTS-11T-6.5-110
115 mm	32 mm	EGTS-11S-6.5-115	EGTS-11T-6.5-115
120 mm	32 mm	EGTS-11S-6.5-120	EGTS-11T-6.5-120
125 mm	32 mm	EGTS-11S-6.5-125	EGTS-11T-6.5-125
130 mm	32 mm	EGTS-11S-6.5-130	EGTS-11T-6.5-130

Cannulated Cancellous $\varnothing 7.3$ -16mm Thread



Core	$\varnothing 4.8$ mm
Outer Diameter	$\varnothing 7.3$ mm
Head Diameter	$\varnothing 9.5$ mm
Head Length	5.6 mm
Drill Bit for thread Hole	5.0 / 2.6
Guide wire	$\varnothing 2.5$ mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	5.0 mm
Tip	Self Tapping / Self Drilling

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
40 mm	16 mm	EGTS-10S-7.3-040	EGTS-10T-7.3-040
45 mm	16 mm	EGTS-10S-7.3-045	EGTS-10T-7.3-045
50 mm	16 mm	EGTS-10S-7.3-050	EGTS-10T-7.3-050
55 mm	16 mm	EGTS-10S-7.3-055	EGTS-10T-7.3-055
60 mm	16 mm	EGTS-10S-7.3-060	EGTS-10T-7.3-060
65 mm	16 mm	EGTS-10S-7.3-065	EGTS-10T-7.3-065
70 mm	16 mm	EGTS-10S-7.3-070	EGTS-10T-7.3-070
75 mm	16 mm	EGTS-10S-7.3-075	EGTS-10T-7.3-075
80 mm	16 mm	EGTS-10S-7.3-080	EGTS-10T-7.3-080
85 mm	16 mm	EGTS-10S-7.3-085	EGTS-10T-7.3-085
90 mm	16 mm	EGTS-10S-7.3-090	EGTS-10T-7.3-090
95 mm	16 mm	EGTS-10S-7.3-095	EGTS-10T-7.3-095
100 mm	16 mm	EGTS-10S-7.3-100	EGTS-10T-7.3-100
105 mm	16 mm	EGTS-10S-7.3-105	EGTS-10T-7.3-105
110 mm	16 mm	EGTS-10S-7.3-110	EGTS-10T-7.3-110
115 mm	16 mm	EGTS-10S-7.3-115	EGTS-10T-7.3-115
120 mm	16 mm	EGTS-10S-7.3-120	EGTS-10T-7.3-120
125 mm	16 mm	EGTS-10S-7.3-125	EGTS-10T-7.3-125
130 mm	16 mm	EGTS-10S-7.3-130	EGTS-10T-7.3-130

Cannulated Cancellous $\varnothing 7.3$ -32mm Thread



Core	$\varnothing 4.8$ mm
Outer Diameter	$\varnothing 7.3$ mm
Head Diameter	$\varnothing 9.5$ mm
Head Length	5.6 mm
Drill Bit for thread Hole	5.0 / 2.6
Guide wire	$\varnothing 2.5$ mm
Screw driver	Cannulated Hexagonal
Hexagonal Size	5.0 mm
Tip	Self Tapping / Self Drilling

Screw Length	Thread Length	S.S. Ref.	Tit. Ref
40 mm	32 mm	EGTS-11S-7.3-040	EGTS-11T-7.3-040
45 mm	32 mm	EGTS-11S-7.3-045	EGTS-11T-7.3-045
50 mm	32 mm	EGTS-11S-7.3-050	EGTS-11T-7.3-050
55 mm	32 mm	EGTS-11S-7.3-055	EGTS-11T-7.3-055
60 mm	32 mm	EGTS-11S-7.3-060	EGTS-11T-7.3-060
65 mm	32 mm	EGTS-11S-7.3-065	EGTS-11T-7.3-065
70 mm	32 mm	EGTS-11S-7.3-070	EGTS-11T-7.3-070
75 mm	32 mm	EGTS-11S-7.3-075	EGTS-11T-7.3-075
80 mm	32 mm	EGTS-11S-7.3-080	EGTS-11T-7.3-080
85 mm	32 mm	EGTS-11S-7.3-085	EGTS-11T-7.3-085
90 mm	32 mm	EGTS-11S-7.3-090	EGTS-11T-7.3-090
95 mm	32 mm	EGTS-11S-7.3-095	EGTS-11T-7.3-095
100 mm	32 mm	EGTS-11S-7.3-100	EGTS-11T-7.3-100
105 mm	32 mm	EGTS-11S-7.3-105	EGTS-11T-7.3-105
110 mm	32 mm	EGTS-11S-7.3-110	EGTS-11T-7.3-110
115 mm	32 mm	EGTS-11S-7.3-115	EGTS-11T-7.3-115
120 mm	32 mm	EGTS-11S-7.3-120	EGTS-11T-7.3-120
125 mm	32 mm	EGTS-11S-7.3-125	EGTS-11T-7.3-125
130 mm	32 mm	EGTS-11S-7.3-130	EGTS-11T-7.3-130

Compression Screws (Herbert)

Design Features

- Hexagonal head for strong and better alignment with screw driver
- Cannulated screw for minimal invasive technique and better trajectory
- Partially threaded for easier compression
- The screw is self-tapping self drilling
- Head with self tapping flutes to facilitate counter sinking of the screw
- Two different pitch of head and shaft threads for easy and controlled compression of the fracture gap
- Reverse cutting distal thread

Material

- Tit

Sizes

- Diameter 2.7mm
 - Length Starting from 10mm to 30 mm with 2mm increments
- Diameter 3.5mm
 - Length Starting from 10mm to 30 mm with 2mm increments
 - Length Starting from 30mm to 50 mm with 5mm increments
- Diameter 5.0mm
 - Length Starting from 16mm to 40 mm with 2mm increments
 - 40mm to 80 mm with 5mm increments
- Diameter 6.5mm
 - Length Starting from 25mm to 120 mm with 5mm increments

Compression Screw $\varnothing 2.7 / \varnothing 3.2$



Screw driver	Hexagonal
Hexagonal Size	1.5 mm
Tip	Self Tapping
Guide wire	

Screw Length	Shaft Length	Thread Length	Tit. Ref
10 mm	2 mm	4 mm	EGTS-12T-2.7-010
12 mm	2 mm	5 mm	EGTS-12T-2.7-012
14 mm	4 mm	5 mm	EGTS-12T-2.7-014
16 mm	4 mm	6 mm	EGTS-12T-2.7-016
18 mm	6 mm	6 mm	EGTS-12T-2.7-018
20 mm	6 mm	7 mm	EGTS-12T-2.7-020
22 mm	8 mm	7 mm	EGTS-12T-2.7-022
24 mm	10 mm	7 mm	EGTS-12T-2.7-024
26 mm	11 mm	8 mm	EGTS-12T-2.7-026
28 mm	12 mm	8 mm	EGTS-12T-2.7-028
30 mm	14 mm	8 mm	EGTS-12T-2.7-030

Compression Screw $\varnothing 3.5 / \varnothing 4.5$



Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping
Guide wire	1 mm / 220 mm

Screw Length	Shaft Length	Thread Length	Tit. Ref
10 mm	2 mm	4 mm	EGTS-12T-3.5-010
12 mm	2 mm	5 mm	EGTS-12T-3.5-012
14 mm	4 mm	5 mm	EGTS-12T-3.5-014
16 mm	5 mm	6 mm	EGTS-12T-3.5-016
18 mm	7 mm	6 mm	EGTS-12T-3.5-018
20 mm	7 mm	7 mm	EGTS-12T-3.5-020
22 mm	9 mm	7 mm	EGTS-12T-3.5-022
24 mm	10 mm	7 mm	EGTS-12T-3.5-024
26 mm	11 mm	8 mm	EGTS-12T-3.5-026
28 mm	12 mm	8 mm	EGTS-12T-3.5-028
30 mm	14 mm	8 mm	EGTS-12T-3.5-030
35 mm	17 mm	9 mm	EGTS-12T-3.5-035
40 mm	18 mm	11 mm	EGTS-12T-3.5-040
45 mm	21 mm	12 mm	EGTS-12T-3.5-045
50 mm	26 mm	12 mm	EGTS-12T-3.5-050

Compression Screw $\varnothing 5.0$ / $\varnothing 6.5$



Screw driver	Cannulated Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping
Guidewire	2.2 / 220

Screw Length	Shaft Length	Thread Length	Tit. Ref
16 mm	5 mm	6 mm	EGTS-12T-5.0-016
18 mm	7 mm	6 mm	EGTS-12T-5.0-018
20 mm	7 mm	7 mm	EGTS-12T-5.0-020
22 mm	9 mm	7 mm	EGTS-12T-5.0-022
24 mm	10 mm	7 mm	EGTS-12T-5.0-024
26 mm	11 mm	8 mm	EGTS-12T-5.0-026
28 mm	12 mm	8 mm	EGTS-12T-5.0-028
30 mm	14 mm	8 mm	EGTS-12T-5.0-030
32 mm	12 mm	8 mm	EGTS-12T-5.0-032
34 mm	14 mm	8 mm	EGTS-12T-5.0-034
36 mm	16 mm	10 mm	EGTS-12T-5.0-036
38 mm	18 mm	10 mm	EGTS-12T-5.0-038
40 mm	18 mm	11 mm	EGTS-12T-5.0-040
45 mm	21 mm	12 mm	EGTS-12T-5.0-045
50 mm	26 mm	12 mm	EGTS-12T-5.0-050
55 mm	29 mm	13 mm	EGTS-12T-5.0-055
60 mm	34 mm	13 mm	EGTS-12T-5.0-060
65 mm	37 mm	14 mm	EGTS-12T-5.0-065
70 mm	42 mm	14 mm	EGTS-12T-5.0-070
75 mm	45 mm	15 mm	EGTS-12T-5.0-075
80 mm	50 mm	15 mm	EGTS-12T-5.0-080

Compression Screw 6.5 / $\phi 7.3$

Screw driver	Cannulated Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping
Guidewire	1.0 / 220



Screw Length	Shaft Length	Thread Length	Tit. Ref
25 mm	9 mm	8 mm	EGTS-12T-6.5-025
30 mm	14 mm	8 mm	EGTS-12T-6.5-030
35 mm	17 mm	9 mm	EGTS-12T-6.5-035
40 mm	18 mm	11 mm	EGTS-12T-6.5-040
45 mm	21 mm	12 mm	EGTS-12T-6.5-045
50 mm	26 mm	12 mm	EGTS-12T-6.5-050
55 mm	29 mm	13 mm	EGTS-12T-6.5-055
60 mm	34 mm	13 mm	EGTS-12T-6.5-060
65 mm	37 mm	14 mm	EGTS-12T-6.5-065
70 mm	42 mm	14 mm	EGTS-12T-6.5-070
75 mm	45 mm	15 mm	EGTS-12T-6.5-075
80 mm	50 mm	15 mm	EGTS-12T-6.5-080
85 mm	53 mm	16 mm	EGTS-12T-6.5-085
90 mm	58 mm	16 mm	EGTS-12T-6.5-090
95 mm	63 mm	16 mm	EGTS-12T-6.5-095
100 mm	68 mm	16 mm	EGTS-12T-6.5-100
105 mm	73 mm	16 mm	EGTS-12T-6.5-105
110 mm	78 mm	16 mm	EGTS-12T-6.5-110
115 mm	83 mm	16 mm	EGTS-12T-6.5-115
120 mm	88 mm	16 mm	EGTS-12T-6.5-120

Plate Locking Screws

Design Features

- Hexagonal head for strong and better alignment with screw driver
- Conical screw head for better alignment of the locking screw in the threaded plate hole
- Large core diameter for better bending and shear strength
- Fully threaded screw
- The screw is self-tapping

Material

- SS 316LVM
- Tit

Sizes

- Diameter 3.5mm
 - Length Starting from 10mm to 30 mm with 2mm increments
30mm to 70 mm with 5mm increments
- Diameter 5.0mm
 - Length Starting from 20mm to 60 mm with 2mm increments
60mm to 120 mm with 5mm increments
- Diameter 6.5mm
 - Length Starting from 25mm to 120 mm with 5mm increments

3.5mm Locking Screw



Core	ø 3.0 mm
Outer Diameter	ø 3.5 mm
Head Diameter	ø 8.0 mm
Head Length	5.0 mm
Drill Bit for thread Hole	ø 3.2 mm
Screw driver	Hexagonal
Hexagonal Size	2.5 mm
Tip	Self Tapping

Length	S.S. Ref.	Tit. Ref
10 mm	EGTS-13S-3.5-010	EGTS-13T-3.5-010
12 mm	EGTS-13S-3.5-012	EGTS-13T-3.5-012
14 mm	EGTS-13S-3.5-014	EGTS-13T-3.5-014
16 mm	EGTS-13S-3.5-016	EGTS-13T-3.5-016
18 mm	EGTS-13S-3.5-018	EGTS-13T-3.5-018
20 mm	EGTS-13S-3.5-020	EGTS-13T-3.5-020
22 mm	EGTS-13S-3.5-022	EGTS-13T-3.5-022
24 mm	EGTS-13S-3.5-024	EGTS-13T-3.5-024
26 mm	EGTS-13S-3.5-026	EGTS-13T-3.5-026
28 mm	EGTS-13S-3.5-028	EGTS-13T-3.5-028
30 mm	EGTS-13S-3.5-030	EGTS-13T-3.5-030
35 mm	EGTS-13S-3.5-035	EGTS-13T-3.5-035
40 mm	EGTS-13S-3.5-040	EGTS-13T-3.5-040
45 mm	EGTS-13S-3.5-045	EGTS-13T-3.5-045
50 mm	EGTS-13S-3.5-050	EGTS-13T-3.5-050
55 mm	EGTS-13S-3.5-055	EGTS-13T-3.5-055
60 mm	EGTS-13S-3.5-060	EGTS-13T-3.5-060
65 mm	EGTS-13S-3.5-065	EGTS-13T-3.5-065
70 mm	EGTS-13S-3.5-070	EGTS-13T-3.5-070

5.0mm Locking Screw



Core	ø 4.0 mm
Outer Diameter	ø 5.0 mm
Head Diameter	Ø 6.5 mm
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Length	S.S. Ref.	Tit. Ref
10 mm	EGTS-13S-5.0-010	EGTS-13T-5.0-010
12 mm	EGTS-13S-5.0-012	EGTS-13T-5.0-012
14 mm	EGTS-13S-5.0-014	EGTS-13T-5.0-014
16 mm	EGTS-13S-5.0-016	EGTS-13T-5.0-016
18 mm	EGTS-13S-5.0-018	EGTS-13T-5.0-018
20 mm	EGTS-13S-5.0-020	EGTS-13T-5.0-020
22 mm	EGTS-13S-5.0-022	EGTS-13T-5.0-022
24 mm	EGTS-13S-5.0-024	EGTS-13T-5.0-024
26 mm	EGTS-13S-5.0-026	EGTS-13T-5.0-026
28 mm	EGTS-13S-5.0-028	EGTS-13T-5.0-028
30 mm	EGTS-13S-5.0-030	EGTS-13T-5.0-030
32 mm	EGTS-13S-5.0-032	EGTS-13T-5.0-032
34 mm	EGTS-13S-5.0-034	EGTS-13T-5.0-034
36 mm	EGTS-13S-5.0-036	EGTS-13T-5.0-036
38 mm	EGTS-13S-5.0-038	EGTS-13T-5.0-038
40 mm	EGTS-13S-5.0-040	EGTS-13T-5.0-040
42 mm	EGTS-13S-5.0-042	EGTS-13T-5.0-042
44 mm	EGTS-13S-5.0-044	EGTS-13T-5.0-044
46 mm	EGTS-13S-5.0-046	EGTS-13T-5.0-046
48 mm	EGTS-13S-5.0-048	EGTS-13T-5.0-048
50 mm	EGTS-13S-5.0-050	EGTS-13T-5.0-050
52 mm	EGTS-13S-5.0-052	EGTS-13T-5.0-052
54 mm	EGTS-13S-5.0-054	EGTS-13T-5.0-054
56 mm	EGTS-13S-5.0-056	EGTS-13T-5.0-056

Screw Length	S.S. Ref.	Tit. Ref
58 mm	EGTS-13S-5.0-058	EGTS-13T-5.0-058
60 mm	EGTS-13S-5.0-060	EGTS-13T-5.0-060
62 mm	EGTS-13S-5.0-062	EGTS-13T-5.0-062
64 mm	EGTS-13S-5.0-064	EGTS-13T-5.0-064
66 mm	EGTS-13S-5.0-066	EGTS-13T-5.0-066
68 mm	EGTS-13S-5.0-068	EGTS-13T-5.0-068
70 mm	EGTS-13S-5.0-070	EGTS-13T-5.0-070
72 mm	EGTS-13S-5.0-072	EGTS-13T-5.0-072
74 mm	EGTS-13S-5.0-074	EGTS-13T-5.0-074
76 mm	EGTS-13S-5.0-076	EGTS-13T-5.0-076
78 mm	EGTS-13S-5.0-078	EGTS-13T-5.0-078
80 mm	EGTS-13S-5.0-080	EGTS-13T-5.0-080
85 mm	EGTS-13S-5.0-085	EGTS-13T-5.0-085
90 mm	EGTS-13S-5.0-090	EGTS-13T-5.0-090
95 mm	EGTS-13S-5.0-095	EGTS-13T-5.0-095
100 mm	EGTS-13S-5.0-100	EGTS-13T-5.0-100
105 mm	EGTS-13S-5.0-105	EGTS-13T-5.0-105
110 mm	EGTS-13S-5.0-110	EGTS-13T-5.0-110
115 mm	EGTS-13S-5.0-115	EGTS-13T-5.0-115
120 mm	EGTS-13S-5.0-120	EGTS-13T-5.0-120

6.5mm Locking Screw Cancellous fully threaded



Core	∅ 4.0 mm
Outer Diameter	∅ 6.5 mm
Head Diameter	∅ 8.0 mm
Screw driver	Hexagonal
Hexagonal Size	3.5 mm
Tip	Self Tapping

Length	S.S. Ref.	Tit. Ref
25 mm	EGTS-13S-6.5-025	EGTS-13T-6.5-025
30 mm	EGTS-13S-6.5-030	EGTS-13T-6.5-030
35 mm	EGTS-13S-6.5-035	EGTS-13T-6.5-035
40 mm	EGTS-13S-6.5-040	EGTS-13T-6.5-040
45 mm	EGTS-13S-6.5-045	EGTS-13T-6.5-045
50 mm	EGTS-13S-6.5-050	EGTS-13T-6.5-050
55 mm	EGTS-13S-6.5-055	EGTS-13T-6.5-055
60 mm	EGTS-13S-6.5-060	EGTS-13T-6.5-060
65 mm	EGTS-13S-6.5-065	EGTS-13T-6.5-065
70 mm	EGTS-13S-6.5-070	EGTS-13T-6.5-070
75 mm	EGTS-13S-6.5-075	EGTS-13T-6.5-075
80 mm	EGTS-13S-6.5-080	EGTS-13T-6.5-080
85 mm	EGTS-13S-6.5-085	EGTS-13T-6.5-085
90 mm	EGTS-13S-6.5-090	EGTS-13T-6.5-090
95 mm	EGTS-13S-6.5-095	EGTS-13T-6.5-095
100 mm	EGTS-13S-6.5-100	EGTS-13T-6.5-100
105 mm	EGTS-13S-6.5-105	EGTS-13T-6.5-105
110 mm	EGTS-13S-6.5-110	EGTS-13T-6.5-110
115 mm	EGTS-13S-6.5-115	EGTS-13T-6.5-115
120 mm	EGTS-13S-6.5-120	EGTS-13T-6.5-120



E-mail:info@egifix.com
website:www.egifix.com